

I.V. House TLC® Elbow Splint with Straps

IV Therapy Elbow Splint / Directions

Single Use Only – Fits either arm.

The TLC Elbow Splint can be applied to the arm before or after inserting IV and connecting to extension set or IV tubing.

1.

Unroll the TLC Splint straps and identify the smooth concave surface of the splint. This concave surface will be placed against the patient's arm.

2.

The "forearm/upper arm" label indicates proper placement. Remove the label before applying the splint.

3.

Place the splint under the patient's arm, with the narrow end of the splint toward the hand and the wide end above the elbow. The point of the elbow should align with the center opening. **We recommend a thumbs-up** application whenever possible.

4.

Wrap the first strap around the forearm. Thread the Velcro® tab into slot and pull it through. Fold the strap around the outer edge of the splint and back onto itself. Secure with Velcro tab.

5.

Repeat with the other two straps, securing each with the Velcro tab.

6.

Note: Slip one finger beneath each strap to ensure the fabric is not so tight that it impedes circulation or infusion of IV fluids.

7. Optional Strap Configuration

If the elbow strap obscures the catheter insertion site, consider a crisscross strap configuration instead. Put the elbow strap through the upper arm slot and the upper arm strap through the elbow slot.

XS 959XS-Ultra TLC® Elbow Splint
Fits most newborns ages 0–4 months,
or weighing 2.5–7 kg
Small splint, shorter straps

S 959S-Ultra TLC® Elbow Splint
Fits most infants ages 4–12 months,
or weighing 6–12 kg
Small splint, longer straps

M 959M-W-Ultra TLC® Elbow Splint
Fits most toddlers ages 12–36 months,
or weighing 9–15 kg
Large splint, shorter straps

L 959L-W-Ultra TLC® Elbow Splint
Fits most children ages 3–10 years,
or weighing 12–41 kg
Large splint, longer straps

Touch / Look / Compare



I.V. HOUSE®
Protection Over and Above®

REV 2025-06-27

I.V. House, Inc. / St. Louis, MO 63146 USA

314-453-9200 / 800-530-0400 / FAX 314-453-9576 / IVinfo@ivhouse.com / www.ivhouse.com